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2.3 Gantry: Removing the rear cover

2.4 Mobile-CT kit (option): Inspections

SI Transportation lock inspected for damage (6M)	<input type="radio"/>	<input type="radio"/>	<input checked="" type="radio"/>
SI Patient table inspected for damage and loose parts (6M)	<input type="radio"/>	<input type="radio"/>	<input checked="" type="radio"/>
SI Gantry inspected for damage, loose parts and screw connections (6M)	<input type="radio"/>	<input type="radio"/>	<input checked="" type="radio"/>

2.5 PHS1A: Cover removal & checking the horizontal drive toothed belt tension

PM Horizontal drive toothed belt tension checked (6M) ==> NOTE: This check should be performed only once during the first maintenance after system/PHS installation or after horizontal drive or toothed belt replacement.	<input type="radio"/>	<input type="radio"/>	<input checked="" type="radio"/>
Value: N.A.			

2.6 Gantry: Cleaning the slip ring

PM Slip ring assembly cleaned (6M)	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>
------------------------------------	----------------------------------	-----------------------	-----------------------

2.7 Gantry: Cleaning the Power Brush Assembly

PM Power brush assembly cleaned (12M)	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>
---------------------------------------	----------------------------------	-----------------------	-----------------------

2.8 Gantry: Checking/replacing the Power Carbon Brushes

PM Power carbon brushes checked (6M)	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>
PM Power carbon brushes replaced (6M)	<input type="radio"/>	<input type="radio"/>	<input checked="" type="radio"/>

2.9 Gantry: Cleaning the Signal Data Brush Block

PM Signal data brush block cleaned (6M)	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>
---	----------------------------------	-----------------------	-----------------------

2.10 Gantry: Checking/replacing the signal data brushes

PM Signal data brushes checked (6M)	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>
PM Signal data brushes replaced (6M)	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>

2.11 Gantry: Checking/replacing the Gantry Air Filters

PM Gantry air filters checked (6M)	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>
PM Gantry air filters replaced (6M)	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>

Completed

2.12 IMS: Cleaning the ICS and IRS

PM ICS and IRS cleaned (6M)

 2.13 IMS: Replacing the IRS air filter

PM IRS air filter replaced (6M)

 2.14 System: Installing the Covers

SIM Covers and ground wire cables installed (6M)

 2.15 Switching on the System**2.16 System: Cleaning the system**

PM System cleaned (6M)

 2.17 System: Quality assurance test

Q Quality assurance test completed (6M)

 2.18 System: Checking the system functionality

Q Functionality of the CT system checked (6M)

 2.19 Maintenance Documentation**2.20 End of Partial Maintenance****3 Full Maintenance (every 12 months)****3.1 Switching on the System****3.2 Gantry: Tilting the Gantry**PM Gantry tilt over the full \pm tilt range executed (12M) **3.3 Switching off the System****3.4 Gantry: Removing the covers****3.5 PHS1A: Removing the covers****3.6 System: Visual inspection according to DIN EN 62353/IEC 62353**

SI Gantry: Visually inspected (12M)

SI Patient table: Visually inspected (12M)

SI System cabinets: Visually inspected (12M)

SI Balancing weights: Visually inspected (12M)

 Completed

SI	System cables & cable routing: Visually inspected (12M)	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>
SI	Accessories: Visually inspected (12M)	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>
SI	Phantoms & phantom holder: Visually inspected (12M)	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>
SI	Care Vision CT (optional): Visually inspected (12M)	<input type="radio"/>	<input type="radio"/>	<input checked="" type="radio"/>
SI	Care Vision CT cables & cable routing (optional): Visually inspected (12M)	<input type="radio"/>	<input type="radio"/>	<input checked="" type="radio"/>

3.7 PHS: Measuring the leakage current

SI	Leakage current (tabletop plate to GND terminal) < 0.1 mA [AC] (12M) Value: 0.00001	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>
SI	Leakage current (tabletop plate to GND terminal) < 0.01 mA [DC] (12M) Value: 0.00005	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>
SI	Leakage current (ECG leads (option) to GND terminal) < 0.1 mA [AC] (12M) Value: N.A.	<input type="radio"/>	<input type="radio"/>	<input checked="" type="radio"/>
SI	Leakage current (ECG leads (option) to GND terminal) < 0.01 mA [DC] (12M) Value: N.A.	<input type="radio"/>	<input type="radio"/>	<input checked="" type="radio"/>
SI	Leakage current (ECG leads (option) to tabletop plate) < 0.1 mA [AC] (12M) Value: N.A.	<input type="radio"/>	<input type="radio"/>	<input checked="" type="radio"/>
SI	Leakage current (ECG leads (option) to tabletop plate) < 0.01 mA [DC] (12M) Value: N.A.	<input type="radio"/>	<input type="radio"/>	<input checked="" type="radio"/>

3.8 Mobile Kit (option): Inspections

SI	Transportation lock inspected for damage (6M)	<input type="radio"/>	<input type="radio"/>	<input checked="" type="radio"/>
SI	Patient table inspected for damage and loose parts (6M)	<input type="radio"/>	<input type="radio"/>	<input checked="" type="radio"/>
SI	Gantry inspected for damage, loose parts and screw connections (6M)	<input type="radio"/>	<input type="radio"/>	<input checked="" type="radio"/>

3.9 System: Inspecting the labels

SI	System labels inspected (12M)	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>
----	-------------------------------	----------------------------------	-----------------------	-----------------------

3.10 PHS1A: Lubricating the vertical Drive Spindle

PM	Vertical drive spindle lubricated (12M)	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>
----	---	----------------------------------	-----------------------	-----------------------

3.11 PHS1A: Lubricating the Guide Rails of the Table

PM	Guide rails lubricated (12M)	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>
----	------------------------------	----------------------------------	-----------------------	-----------------------

3.12 PHS1A: Checking the Horizontal Drive Toothed Belt Tension

PM	Horizontal drive toothed belt tension check (6M) - A check of the drive belt tension should be performed only once at the first maintenance after system/ PHS installation or after horizontal drive or toothed belt replacement. Value: N.A.	<input type="radio"/>	<input type="radio"/>	<input checked="" type="radio"/>
----	--	-----------------------	-----------------------	----------------------------------

Completed

3.13 PHS1A: Checking the anchoring of the PHSPM Anchoring checked (12M) **3.14 MCS (option): Checking the monitor ceiling system**SI MCS visually checked for varnish damage (12M) SI MCS visually checked for cracks (12M) PM MCS mechanical movements checked (12M) PM MCS load balance checked (12M) **3.15 Gantry: Checking the rotation drive belt and belt tension**PM Belt checked (12M) PM Belt tension checked (12M) **3.16 Gantry: Checking the Rotation Brake**PM Rotation brake checked (12M) **3.17 Gantry: Checking the Cover Hooks**PM Elasticity of the hook holders on the front and rear gantry side checked (12M) **3.18 Gantry: Cleaning the Air Outlets on the Gantry Top**PM Gantry air outlets cleaned (12M) **3.19 Gantry: Checking the Gantry Cooling Fans**PM Gantry cooling fans checked (12M) **3.20 Gantry: Cleaning the slip ring**PM Slip ring assembly cleaned (12M) **3.21 Gantry: Cleaning the Power Brush Assembly**PM Power brush assembly cleaned (6M) **3.22 Gantry: Checking/replacing the Power Carbon Brushes**PM Power carbon brushes checked (6M) PM Power carbon brushes replaced (6M) **3.23 Gantry: Cleaning the Signal Data Brush Assembly**PM Signal data brush assembly cleaned (6M) **Completed**

3.24 Gantry: Checking/replacing the signal data brushes

PM Signal data brushes checked (6M)

PM Signal data brushes replaced (6M)

 3.25 Gantry: Checking/replacing the Gantry Air Filters

PM Gantry air filters checked (6M)

PM Gantry air filters replaced (6M)

 3.26 Gantry: Lubricating the Main Bearing

PM Main bearing lubricated (12M)

 3.27 Gantry: Checking/replacing the MAS Air Filter

PM MAS air filter checked (12M)

PM MAS air filter replaced or washed (12M)

 3.28 Gantry: Checking/Replacing the DMS Air Filter

PM DMS air filter checked (12M) - only for 64/128 slice systems

PM DMS air filter replaced (12M) - only for 64/128 slice systems

 3.29 Gantry: Cleaning the Tube Cooler

PM Tube cooler cleaned (12M)

 3.30 IMS: Cleaning the ICS and IRS

PM ICS and IRS cleaned (6M)

 3.31 IMS: Replacing the Air Filter for IRS G1/G2

PM IRS air filter replaced (6M)

 3.32 System: Measuring the protective conductor resistance**3.32.2 SOMATOM Perspective**SIE LCB: ground wire resistance $\leq 100 \text{ m}\Omega$ (12M)

Value: 10

SIE Gantry: ground wire resistance $\leq 300 \text{ m}\Omega$ (12M)

Value: 40

SIE Patient table: ground wire resistance $\leq 300 \text{ m}\Omega$ (12M)

Value: 60

SIE ICS Tower PC: ground wire resistance $\leq 300 \text{ m}\Omega$ (12M)

Value: 60

Completed

		OK	not OK	n.a.
SIE ICS monitor: ground wire resistance $\leq 300 \text{ m}\Omega$ (12M)		<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>
Value:	80			
SIE UPS for ICS: ground wire resistance $\leq 300 \text{ m}\Omega$ (12M)		<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>
Value:	60			
SIE IRS Tower PC: ground wire resistance $\leq 300 \text{ m}\Omega$ (12M)		<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>
Value:	70			
SIE IES Tower PC (option): ground wire resistance $\leq 300 \text{ m}\Omega$ (12M)		<input type="radio"/>	<input type="radio"/>	<input checked="" type="radio"/>
Value:	N.A.			
SIE IES monitor (option): ground wire resistance $\leq 300 \text{ m}\Omega$ (12M)		<input type="radio"/>	<input type="radio"/>	<input checked="" type="radio"/>
Value:	N.A.			
SIE UPS for IES (option): ground wire resistance $\leq 100 \text{ m}\Omega$ (12M)		<input type="radio"/>	<input type="radio"/>	<input checked="" type="radio"/>
Value:	N.A.			
SIE Ceiling support (in scan room, option): ground wire resistance $\leq 300 \text{ m}\Omega$ (12M)		<input type="radio"/>	<input type="radio"/>	<input checked="" type="radio"/>
Value:	N.A.			
SIE Monitor trolley (in scan room, option): ground wire resistance $\leq 300 \text{ m}\Omega$ (12M)		<input type="radio"/>	<input type="radio"/>	<input checked="" type="radio"/>
Value:	N.A.			
SIE Monitor (in scan room, option): ground wire resistance $\leq 300 \text{ m}\Omega$ (12M)		<input type="radio"/>	<input type="radio"/>	<input checked="" type="radio"/>
Value:	N.A.			

3.33 Switching on the System

3.34 System: Function Tests

SI Emergency STOP circuit: Function checked (12M)	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>
SI Radiation cut-off: Function checked (12M)	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>
SI Radiation indicators: Function checked (12M)	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>
SI Door contact switch: Function checked (12M)	<input type="radio"/>	<input type="radio"/>	<input checked="" type="radio"/>
SI Patient positioning aids: Function checked (12M)	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>
SI Line connection box: Function checked (12M)	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>

3.35 Gantry: Checking the DMS Fans

PM DMS Fans checked (12M)	<input type="radio"/>	<input type="radio"/>	<input checked="" type="radio"/>
---------------------------	-----------------------	-----------------------	----------------------------------

3.36 System: Reinstalling the covers

SIE Covers and protective conductor cables installed (6M)	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>
---	----------------------------------	-----------------------	-----------------------

3.37 System: Cleaning the system

PM System cleaned (6M)	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>
------------------------	----------------------------------	-----------------------	-----------------------

3.38 System: Quality assurance test

Q Quality assurance test completed (6M)	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>
---	----------------------------------	-----------------------	-----------------------

Completed

3.39 System: Checking the System Functionality

Q Functionality of the CT system checked (6M)

3.40 Maintenance Documentation

3.41 End of full maintenance

4 Preventive Parts Replacement

4.1 IVM (optional): Replacing the Batteries in the IVP

PM IVP batteries replaced (24 M)

Startup Date:

N.A. · N.A. · N.A.

Date of Last Replacement:

DD · MMM · YYYY

4.2 IMS: Replacing the UPS battery pack

PM Battery pack of UPS replaced (30M)

Startup Date:

18 · Jan · 2023

Date of Last Replacement:

DD · MMM · YYYY

4.3 IES: Replacing the complete UPS

PM Complete UPS replaced (24M)

Startup Date:

N.A. · N.A. · N.A.

Date of Last Replacement:

DD · MMM · YYYY

4.4 PHS1A: Replacing the Gas Spring

PM Gas spring replaced (42M)

Startup Date:

19 · Jul · 2022

Date of Last Replacement:

DD · MMM · YYYY

4.5 End of Preventive Parts Replacement

Completed